LOGS

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The image shows a chemical structure editor interface. At the top, there is a toolbar with icons for copy, paste, save, undo, redo, clear, and delete. On the left side, there is a vertical toolbar with icons for mouse selection, drawing a bond, drawing a line, a carbon atom (C), a ring structure, and a plus/minus sign. The main area is a large, empty white space for drawing the chemical structure. At the bottom, there is another toolbar with icons for zoom in, zoom out, pan, and other navigation functions.

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


The screenshot displays a software interface for scientific data management. On the left, a list of items is shown, with item 8 selected. The interface includes a toolbar with buttons for Submit, Cancel, Edit, Delete, Add, Copy, and a numeric input field. On the right, a 'Title' dialog box is open, featuring a text input field and two buttons: 'Set Title' and 'Set & Copy Title'.

Item ID	Status	Item Name	Quantity	Unit
8	Available	JEP-132-5	10	N 1H
9	Available			
10	Available			
11	Available			
12	Available			

Buttons: Submit, Cancel, Edit, Delete, Add, Copy, 1, 1


Dialog Box: Title, Set Title, Set & Copy Title

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
Integration in TopSpin

The screenshot shows the TopSpin software interface with the 'AcqPars' tab selected. The 'Automation' section is highlighted with a red box. The parameters are as follows:

Parameter	Value
LOCNOC	ZH
SOLVENT	DMSO
Automation	
AUNM	au_zg
EXP	1H
Miscellaneous	
GRDPROG	
CHEMSTR	none
User parameters	
USERA1	user
USERA2	user
USERA3	user

```
GETCURDATA
GETPROSOL
XCMD("sendgui xpy getMolFileFromWeb.py")
RGA
ZG
QUIT
```

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Experiment Table

Ho...	Type	Status	Name	No.	Solvent	Experiment	Pri	Par	Title/Orig
32	2	Finished							
	Finished	IDP-PBZ				N 1H	★		IDP-PBZ #202391# [redacted] 1H
	Finished	IDP-PBZ				N 31P	★		IDP-PBZ #202391# [redacted] 31P
33	1	Finished							
	Finished	IDP-PB3				N 1H	★		IDP-PB3 #202391# [redacted]

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




Dataset #1890-2 13C{1H} - MLK-J08B ▶ MLK-J08B

Data Origin: data-server//home/logstransfer/nmrdata/av250-oc/AG_X/nmr/MLK-J08B/11


[Download](#) [Edit](#) [Trash](#)

Entry: NMR - Thu, Jan 28, 2021

[Details](#) [Documents](#) [More data on this sample](#) [File contents](#) [View spectrum](#)

Zoom: [Reset](#) Tools: Track slider Curve stacking Processing
 Select box Auto zoom

[Export](#)
[Settings](#)

- MLK-J08B/11
 - spectra
 - proc 1
 - fid
 - fid.re
 - fid.im
 - LOGS generated
 - spectrum

Acquisition parameters


- Pulse program zgpg30
- Scans 1024
- Dummy scans 4
- Receiver gain 2048
- Spectral width 451.637114206743 ppm
- Spectral width 28409.0909090909 Hz
- Acquisition time 1.153 s
- Size of FID 65536
- Transmitter frequency 62.895239 MHz
- Frequency offset 115.0 ppm
- Temperature 297.0163 K
- Probe 5 mm BBO BB-1H Z-GRD Z3074/
- MAS Spin Rate 4200

Dataset #720-7 HSQCFFP3GPPHWG - 1H, 15N HSQC (with 13C Decoupling) ▶ PB-Ubiquitin

Data Origin: <data-server//home/logstransfer/nmrdata/examples/Research/nmr/PB-Ubiquitin/5>

[Download](#)[Edit](#)[Trash](#)

Entry: NMR - Mon, Jan 18, 2021

[Details](#) [Documents](#) [More data on this sample](#) [File contents](#) [View NMR spectrum \(Bruker\)](#)[Spectrum](#) [FID](#) [fft on f2 dimension](#) [Bruker Acqus File](#) [Bruker Acqu2s File](#) [Bruker Pulseprogram](#) hide negative

contour level

contour factor

 phase rows phase columns

Filter

Acquisition parameters

Pulse program hsqcftp3gpphwg

Scans 4

Dummy scans 16

Receiver gain 451.44

Spectral width 16.0215603155368 ppm, 34.997909136943 ppm

Spectral width 9615.38461538462 Hz, 2128.56534695615 Hz

Acquisition time 106.496 ms, 60.134 ms

Size of FID 2048, 256

Transmitter frequency 600.15 MHz, 60.812671 MHz

Frequency offset 4.69999166874948 ppm, 116.99995877503883 ppm

Temperature 298.0713 K

Probe Z856301_0002 (PA QXI 600S3 H/P-C/N-D-05 XYZ BTO)

MAS Spin Rate 4200

Nucleus 1H, 15N

General information

Measurement duration 20 min 5 s

Dimension 2

LOGS Projects Samples Datasets Documents Search Monitoring More

Dataset #720-7 HSQCFF3GPPHWG - 1H, 15N HSQC
 Data Origin: data-server//home/logstransfer/nmrdata/examples/Research/nmr/PB-Ubiqutin/5

Details Documents More data on this sample File contents

Spectrum FID fft on f2 dimension Bruker Acqus File Bruker Acqus2s File Bruker

Phasing F2 at 127.28 ppm (row 52) Save phase Cancel

0th order correction

1st order correction

12.71

$\delta(1H)$ [ppm]

Pivot point (ppm)

contour level

contour factor


Filter

Acquisition parameters

- Pulse program hsqcff3gpphwg
- Scans 4
- Dummy scans 16
- Receiver gain 451.44
- Spectral width 16.0215603155368 ppm, 34.997909136943 ppm
- Spectral width 9615.38461538462 Hz, 2128.56534695615 Hz
- Acquisition time 106.496 ms, 60.134 ms
- Size of FID 2048, 256
- Transmitter frequency 600.15 MHz, 60.812671 MHz
- Frequency offset 4.69999166874948 ppm, 116.99995877503883 ppm
- Temperature 298.0713 K
- Probe Z856301_0002 (PA QXI 600S3 H/P-C/N-D-05 XYZ BTO)
- MAS Spin Rate 4200
- Nucleus 1H, 15N

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Guidelines for Safeguarding Good Research Practice

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Kodex

<https://doi.org/10.5281/zenodo.6472827>

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„[...] all universities and non-university research institutions will be required to implement the 19 guidelines in a legally binding manner to be eligible to receive DFG funding.“

The four-year transition period for implementing ends 31 July 2023

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Guidline 17 – Archiving: Researchers back up research data and results made publicly available, as well as the central materials on which they are based and the research software used, by adequate means according to the standards of the relevant subject area, and retain them for an appropriate period of time [10 years]. The archiving period begins on the date when the results are made publicly available.

SCIENTIFIC DATA

Amended: Addendum

OPEN

SUBJECT CATEGORIES

- » Research data
- » Publication characteristics

Comment: The FAIR Guiding Principles for scientific data management and stewardship

Wilkinson, M. D. et al. The FAIR Guiding Principles for scientific data management and stewardship. *Sci Data* **3**, 160018 (2016).

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Findable Accessible Interoperable Reusable

To be Findable:

- F1. (meta)data are assigned a globally unique and persistent identifier
- F2. data are described with rich metadata (defined by R1 below)
- F3. metadata clearly and explicitly include the identifier of the data it describes
- F4. (meta)data are registered or indexed in a searchable resource

To be Accessible:

- A1. (meta)data are retrievable by their identifier using a standardized communications protocol
 - A1.1 the protocol is open, free, and universally implementable
 - A1.2 the protocol allows for an authentication and authorization procedure, where necessary
- A2. metadata are accessible, even when the data are no longer available

Wilkinson, M. D. et al. The FAIR Guiding Principles for scientific data management and stewardship. *Sci Data* **3**, 160018 (2016).

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Findable Accessible Interoperable Reusable

To be Interoperable:


- I1. (meta)data use a formal, accessible, shared, and broadly applicable language for knowledge representation.
- I2. (meta)data use vocabularies that follow FAIR principles
- I3. (meta)data include qualified references to other (meta)data

To be Reusable:

- R1. meta(data) are richly described with a plurality of accurate and relevant attributes
 - R1.1. (meta)data are released with a clear and accessible data usage license
 - R1.2. (meta)data are associated with detailed provenance
 - R1.3. (meta)data meet domain-relevant community standards

Wilkinson, M. D. et al. The FAIR Guiding Principles for scientific data management and stewardship. *Sci Data* **3**, 160018 (2016).

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```
+++ Result: +++
```

```
User 'Baker, Charles' measured 42 datasets on 22 samples  
the total Measuring time was 22800 s (6:20:00)  
with a measuring cost of 120 €/h this costs 760.0 €
```

```
User 'Parker, Nancy' measured 20 datasets on 20 samples  
the total Measuring time was 1200 s (0:20:00)  
with a measuring cost of 120 €/h this costs 40.0 €
```

bill_000002.html

30,6 kB

LOGS
Python

bill_000002.pdf

87,5 kB

LOGS

Python

Invoice 000002

Baker, Charles

charles.baker@davis.com

Samples (22): #8 Pepsi , #9 Pepsi Light , #4 Coca Cola Zero Koffeinfrei , #18 Afri Cola , #19 Coca Cola (USA) , #7 Dr. Pepper , #20 Coca Cola (Mexico) , #11 Sinalco Cola , #14 Bizzl Cola , #17 Fritz Cola , #3 Coca Cola Zero , #21 Ubiquitin , #5 Coca Cola Vanille , #48 Sample Name , #13 Korrekt Cola , #2 Coca Cola Light , #12 Bionade Cola , #15 Vita Cola , #16 Red Bull Cola , #1 Coca Cola , #10 Pepsi Max , #6 Fanta Orange

Datasets measured: 42

Measurement time: 6:20:00

Total cost: 760.0€ (cost per hour 120€)

Created at 22/03/2023 08:17:20

bill_000003.html

28,0 kB

LOGS

Python

bill_000003.pdf

84,7 kB

LOGS

Python

Invoice 000003

Campbell, Nicole

nicole.campbell@harris.edu

Samples (1): #43 Azol

Datasets measured: 9

Measurement time: 2 days, 11:48:33

Total cost: 7177.1€ (cost per hour 120€)

Created at 22/03/2023 08:17:20

LOGS

Search...

Contents

When Version Who What

2021-07-16 0.1 Sina Kazemi First version

API Documentation

class LOGS

- LOGS()
- getHumanReadableSize()
- update()
- create()
- delete()
- deleteById()
- sample()
- samples()
- project()
- projects()
- dataset()
- datasets()
- document()
- documents()
- person()

LOGS.LOGS

Purpose: A library to access the LOGS API via Python Python Version: Python 3.8 and above Disclaimer: This script is provided "as is" by Signals GmbH & KO KG to demonstrate the functionality of the LOGS API License: MIT License

When Version Who What

2021-07-16 0.1 Sina Kazemi First version

class LOGS:

Python class to access the LOGS web API

```
LOGS(  
    url: str,  
    apiKey: str,  
    use_internal: bool = False,  
    verbose: bool = False,  
    verify: bool = True  
)
```

Checks the connection to the server on creation

Args: url (str): URL to specific LOGS group (e.g. <https://mylogs/mygroup> or <https://mylogs:80/mygroup/api/0.1>) api_key (str): The API key that grants access to LOGS (you need to generate one in LOGS and copy it) verbose (bool, optional): If set you see some information about the server connection. Defaults to False.

Raises: Exception: URL does not defined or is invalid. Exception: The URL does not define a group. Exception: Server cannot be reached.

```
@classmethod  
def getHumanReadableSize(cls, size, suffix='B'):
```

LOGS

Python

Example application cases

Question for the audience:

Which type of application examples could you think of?

What kind of automated data extraction or manipulation could help you with your daily work?

(e.g. routine tasks, monitoring of properties, automated data manipulation and analysis)

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